

BioXC™ 200GX

TRIGGERED OR CONTINUOUS AIR SAMPLER

The BioXC™ 200GX is an affordable approach to continuous or triggered bio-threat sampling. It is designed for seamless integration with the Cepheid GeneXpert® presumptive identifier. The innovative design of the BioXC eliminates the need to transfer samples from a collection vial or cartridge into the GeneXpert cartridge, significantly reducing risk of inadvertent human exposure and re-aerosolization of concentrated pathogens.

The approach helps speed up the time from detection to confirmation testing and minimize the risk of human contact with the sample. When a GeneXpert cartridge is placed into the BioXC dispensing chamber, the sample is automatically deposited into the cartridge and ready for subsequent analysis. A standard sample vial can be used similarly.

By automatically dispensing the liquid into the cartridge, the BioXC eliminates the time-consuming and costly process of transferring the sample manually with a pipette. No additional fluidics are required.

The BioXC may also be integrated with the FLIR AirSentinel® 1000B to create a trigger-based bio-threat detection system. In that scenario, the AirSentinel 1000B would trigger the BioXC to take a sample only under threatening conditions, rather than continuously sampling, thus reducing the costs of consumables.



FEATURES


- 150 lpm air flow rate
- Continuous autonomous operation
- Optional triggered sample collection
- Recommended for use with the FLIR AirSentinel 1000B detector
- May be integrated with third-party biological identification and detection technology
- Flexible design accommodates unique customer requirements and enables a variety of deployment models
- Optional wireless communications

The BioXC uses the FLIR patented rotating impactor technology. This unique air-sweep approach allows high-volume aerosol samples with low power requirements and reduced risk of reaerosolization relative to wet-walled cyclones.

The BioXC isolates the aerosol sample into a standard vial or the GeneXpert cartridge. The auto-dispense feature offers a safe and simple approach to sample collection. To initiate the testing sequence, remove the cartridge from the dispensing station and place it into the GeneXpert presumptive identifier. If a standard vial is used, the liquid sample is ready for analysis.

SPECIFICATIONS

Height	10 in (25.4 cm)
Width	6.8 in (17.2 cm)
Depth	5.5 in (14 cm)
Weight	5 lbs (2.3 kg)
Power Source	24 VDC
Sample Flow Rate	150 lpm
Particle Size Collection	1 to 10 microns
Collection Efficiency	50% average
Operating Temperature	0°C to 40°C
Storage Temperature	-40°C to 70°C
Interface	RS-232
Maximum Power	33.9 watts
Standby Power	2.4 watts



APPLICATIONS

- Building facility management
- Subways and airports
- Mailrooms and mail distribution facilities
- Outdoor event surveillance
- Hospitals and other healthcare facilities
- Educational institutions
- Military

