

Assure

Rapid Antibody Test



KEY FEATURES & BENEFITS

- Fast results (<15 minutes)
- High Accuracy
- Easy to use
- No equipment required

First FDA Authorized rapid test for use in Point-of-Care with a fingerstick specimen

Rapid Antibody Test Kit

The Assure Rapid antibody test uses a fingerstick blood specimen to provide a rapid result for individuals who want to understand their prior exposure to the COVID-19 virus.

Assure Rapid Antibody Test for Point of Care Use

The Assure COVID-19 IgG/IgM Rapid Test Device is intended for use as an aid in identifying individuals with an adaptive immune response to SARS-CoV-2, indicating recent or prior infection.

The Assure Rapid Test has been Emergency Use Authorized for use with fingerstick whole blood specimens at the Point of Care (POC), i.e., in patient care settings operating under a CLIA Certificate of Waiver, Certificate of Compliance, or Certificate of Accreditation.

Important Information

At this time, it is unknown for how long antibodies persist following infection and if the presence of antibodies confers protective immunity. The Assure COVID-19 IgG/IgM Rapid Test Device should not be used to diagnose acute SARS-CoV-2 infection.

The test is for professional use only.

This test is also authorized for use with fingerstick whole blood specimens only at the Point of Care (POC), i.e., in patient care settings operating under a CLIA Certificate of Waiver, Certificate of Compliance, or Certificate of Accreditation.