

# DosiFFR

Dosimeter Management Software  
For First Responders



## OVERVIEW

DosiFFR is a streamlined Windows-based software application that enables users to manage dosimeter, personnel and exposure data. The single workstation system provides for straightforward log-in/log-out of personnel and the configuration of DMC2000, SOR and DMC3000 electronic self-reading dosimeters.

It has been designed to meet the needs of smaller sized organizations that require a simple and user-friendly solution for issuing dosimetry and managing and collecting data. It is ideal for use in the field and in a limited user environment where quick instrument issue and return is a necessary and significant feature.

DosiFFR is the perfect solution for first responders, hazmat teams, laboratories and industrial users looking for a simple, user-friendly way to configure a dosimeter based on their specific requirements.

Simplified Dosimeter Issue and  
Data Collection System



## KEY FEATURES

- Fast Login/logout
- Simple dosimeter activation
- Instant activation with single press of a button
- No reader required for field activation of the dosimeter
- Dosimeter profile feature for dosimeter configuration
- Data management and report generation - DosiExporter
- Data transfer between several DosiFFR Installations (optional)
- User-friendly interface
- Touch screen support for mobile deployment

# DosiFFR


## FUNCTION

- **Dosimeter preparation, issue and return**
  - Batch configuration process
  - Dosimeter settings (thresholds, parameters, assigned user)
  - Activation of fast-entry mode
  - Simple dosimeter return
  - Batch return process
  - Dose assignment to one or several users
  - Reset and prepare dosimeter for next use
  - Manual dose assignment
- **Users Management**
  - Dosimeter profile management
  - Dosimetry visualization and analysis capabilities
- **Dose extraction and reporting**
  - ✓ Optional centralization capabilities through a dosimetry data extraction and transfer from DosiFFR to DosiServ
  - ✓ Customizable users dosimetry data summary (pdf format)
  - ✓ Multi-units dose and rate display



## SYSTEM CONFIGURATION

- Operating system: Microsoft Windows XP and later
- Embedded database for each installation
- (SQL server - free version)
- DMC2000, SOR, DMC3000 dosimeter compliant
- Optional data extraction and reporting tool; DosiExporter



### RELATED PRODUCTS

Dosimeters:  
DMC 2000, SOR, DMC 3000

Readers:  
LDM 220, LDM 320 D, LDM 320 W