

Bluetooth® Adaptor And Mobile App RadEye/PackEye Accessory



Faster response to alarms without exposing operation

- Free app for your Apple® or Android™ device
- Easily install the Bluetooth adaptor to your Thermo Scientific™ RadEye™ or Thermo Scientific™ PackEye™ without tools
- Upgrade your instrument in minutes with no tools required Real time alarm indications on every connected smartphone
- See data on your phone in real time
- Receive alarms and rate data on your connected smart watch
- See detailed radiological and system information on your mobile device including count rate, dose rate, and alarms.



RadEye




PackEye



Applicable RadEye models are RadEye RadEye SPRD/GN RadEye PRD and RadEye PRD-ER.



Alarm activated on RadEye SPRD App



Natural radiation detected – no threat to public

Download a copy of our apps today.
Available for iOS and Android operating systems.



Ordering Information	
RadEye BT Adaptor	425067095
PackEye BT4 Adaptor	SM168535309

